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
REFERENCES

Atlas RM, Williams JF, Huntington MK. Legionella contamination of dental-unit waters. Applied and Environmental Microbiology. 1995 Apr 1;61(4):1208-13.

Arvand M¹, Hack A. Eur J Microbiol Immunol (Bp). Microbial contamination of dental unit waterlines in dental practices in Hesse, Germany: A cross-sectional study. 2013 Mar;3(1):49-52. doi: 10.1556/EuJMI.3.2013.1.7. Epub 2013 Mar 13.

Castiglia P, Liguori G, Montagna MT, Napoli C, Pasquarella C, Bergomi M, Fabiani L, Monarca S, Petti S, SItI Working Group Hygiene in Dentistry Italian multicenter study on infection hazards during dental practice: control of environmental microbial contamination in public dental surgeries. BMC Public Health. 2008 May 29;8:187. doi: 10.1186/1471-2458-8-187.

Fotos PG, Westfall HN, Snyder IS, Miller RW, Mutchler BM. Prevalence of Legionella-specific IgG and IgM antibody in a dental clinic population. Journal of Dental Research. 1985 Dec 1;64(12):1382-5.

Pankhurst CL, Johnson NW, Woods RG. Microbial contamination of dental unit waterlines: the scientific argument. International Dental Journal. 1998 Aug 1;48(4):359-68.

Ricci ML, Fontana S, Pinci F, Fiumana E, Pedna MF, Farolfi P, Sabattini MA, Scaturro M. Pneumonia associated with a dental unit waterline. Lancet. 2012 Feb 18;379(9816):684. doi: 10.1016/S0140-6736(12)60074-9.

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Von Baum H, Ewig S, Marre R, Suttorp N, Gonschior S, Welte T, Lück C, Competence Group Community-acquired Legionella pneumonia: new insights from the German competence network for community acquired pneumonia. Clinical infectious diseases: an official publication of the Infectious Diseases Society of America. 2008 May 1;46(9):1356-64. doi: 10.1086/586741.

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Unit-Oral Sampling kit

Unit-Oral Sampling Kit for Water Stabilization

INSTRUCTIONS FOR USE

UNIT-ORAL

Sampling Kit for Dental Unit Water

Kit Contents :

1 tube with reference label, 1 protection container, 1 prescription form, 1 postal envelope and instructions for use.

Precautions :

Perform water collection before or 2 days after the last treatment of the dental unit waterlines with antibacterial solution. Sterile gloves should be worn when collecting. Close tightly the tube and the container to avoid water leak.

Instructions for use :

1) Wash the outside of the turbine with antibacterial solution or with antibacterial wipe.

2) Turn on the turbine without using it for at least two minutes to drain the waterlines.

3) Fill the provided tube with water. The tube must absolutely be full (above the line).

4) Close tightly the tube. Check the reference label is stuck on the tube.

5) Put the tube in the protection container to secure the transport of the sample.

6) Send the sample and the fulfilled prescription form to the laboratory with the postal envelope.

INSTRUCTIONS FOR USE

UNIT-ORAL

Sampling Kit for stabilization of Dental Unit Water and DNA/RNA analysis

1) Perform water collection before or 2 days after the last treatment with antibacterial solution. Wash the outside of the turbine with antibacterial solution or with antibacterial wipe. Sterile gloves should be worn when collecting.

2) Turn on the turbine without using it for at least two minutes to drain the waterlines.

3) Unscrew the provided beaker and fill it with 10 ml of water. Screw lid back onto the beaker.

4) Hold beaker on a flat surface and remove the round safety sticker. Take one sample tube (DO NOT OPEN) and firmly insert it in the opening with the cap at the bottom. The sample tube fills automatically. Bring out the sample tube as soon as it is filled.

Repeat the same operation with the second sample tube. Homogenize the contents of the sample tubes and preserve at room temperature until forwarding. After use, the cover in the lid of the beaker must be sealed again using the round safety sticker.

5) Put the tubes in the protection envelope to secure the transport of the sample and send the sample and the fulfilled prescription form to the laboratory with the postal envelope.

Warning statements and precautions:

Beaker: The round safety sticker placed on the lid should be withdrawn only for the transfer of the solution in the sample tube. Then, the sticker must be replaced. **Never** put fingers into the opening because there is a risk of needlestick!

Sample tubes: The tubes contain sodium azide, which is poisonous*. These tubes should not be opened. In case of accidental swallowing, doctor's advice should be sought immediately, with reference to the instructions for use.

*Special precautions for people and environment: R 28 is very poisonous if swallowed. R 32 produces poisonous vapours on contact. R 50/53 is very poisonous for water organisms and can cause long-term damage in waters.